




Sepax Technologies, Inc.

Preparative Size Exclusion Chromatography

SRT[®]-10/10C (10 μ m)

*Fast Protein Purification
For FPLC, HPLC and Prep LC*

*Better Surface Chemistry
for Better Separation*





SRT[®]-10/10C Size Exclusion Columns

FAST Speed, High Resolution and High Loading
Semi-prep and Preparative Size Exclusion Separation
For FPLC, HPLC and Prep LC

General Description

Utilizing proprietary innovative surface technologies, SRT-10/10C SEC phases are made of uniform, hydrophilic, and neutral nanometer thick films chemically bonded on high purity and mechanically stabilized 10 μm silica. The combination of the two different types of SEC coating chemistries, SRT-10 with stand-up monolayer bonded on porous silica and SRT-10C with lay-down monolayer on porous silica (Figure 1), provides a powerful total solution for robust, reproducible and high resolution size based separation of biological molecules in the market.

Our unique bonding chemistry, coupled with the maximized bonding density, allows SRT-10/10C phases to provide high stability and negligible non-specific interactions. SRT-10/10C packings have large pore volumes, resulting in high separation resolution. The columns can tolerate high backpressure, which allows for higher flow rate/faster speed to save the total run time, thus increasing purification efficiency at the laboratory preparative scale. The well-controlled surface chemistry results in excellent lot-to-lot reproducibility.


Typical applications for SRT-10/10C columns include separation and detection of biological molecules like antibodies, antibody drug conjugates, proteins, peptides, oligonucleotides and virus like particles, as well as water-soluble polymers in aqueous buffers.

Featured Characteristics

- Particle size: 10 μm
- Pore size selection: 100, 300, 500, and 1000 \AA
- pH stability 2-8.5
- High capacity and loadability
- Fast speed and high resolution
- High stability over low and high concentration salt
- Excellent lot-to-lot reproducibility
- High protein recovery with intact biological activity
- Chemically bonded stationary phase resulting in negligible non-specific interactions
- Ideal for separation of biological molecules: antibody, proteins, peptides, oligonucleotides, peptides and virus
- Ideal for separation of natural polymers, e.g. polysaccharides, water soluble synthetic polymers, and nanomaterials, e.g. nanoparticles

Stationary Phase Structure

Figure 1. Phase structure difference: a monolayer stands up on the silica surface for SRT-10, and a monolayer lays down on the silica surface for SRT-10C.




Pore size vs. MW exclusion limit

The narrowly dispersed, spherical silica particles of the SRT-10/10C SEC have nominal pore sizes at 100, 300, 500, and 1000 \AA , respectively. Figure 2 and 3 show the separation profiles of protein mixture on the SRT-10 SEC columns with different pore sizes (100, 300, 500, 1000 \AA).


| Phases (10 μm) | Pore Size | Protein MW Exclusion Limit (Da) |
|----------------------------|-------------------|---------------------------------|
| SRT-10/10C SEC-100 | 100 \AA | 30,000 |
| SRT-10/10C SEC-300 | 300 \AA | 1,250,000 |
| SRT-10/10C SEC-500 | 500 \AA | 5,000,000 |
| SRT-10/10C SEC-1000 | 1000 \AA | 7,500,000 |

Figure 2. Comparison of the separation profiles of a protein mixture by SRT-10 SEC columns with different pore sizes: 300, 500 and 1000 \AA , respectively



Columns: SRT-10 SEC (10 μ m, 7.8 x 300 mm)
 Mobile phase: 150 mM phosphate buffer, pH 7.0
 Flow rate: 1.0 mL/min
 Detection: UV 214 nm
 Temperature: Ambient (23 °C)
 Injection volume: 3.0 μ L
 Sample: 1) Thyroglobulin 670 kD; 2) BSA dimer 132 kD; 3) BSA 66 kD; 4) Ribonuclease A 13.7 kD; 5) Uracil 120 D.

Figure 3. The separation profiles of small proteins by SRT-10 SEC 100 Å pore size column



Columns: SRT-10 SEC-100 (10 μ m, 100 Å 21.2 x 300 mm)
 Mobile phase: 150 mM phosphate buffer, pH 7.0
 Flow rate: 7.0 mL/min
 Detection: UV 280 nm
 Temperature: Ambient (23 °C)
 Injection volume: 1.5 mL
 Sample: 3.3 mg/mL BSA 66 kD; 3.3 mg/mL Ribonuclease A 13.7kD; 0.3 mg/mL Uracil 120 D


Fast Purification with High Flow Rate

Sepax SRT-10/10C preparative size SEC columns can run under a faster flow rate than the GE Superdex columns, thus greatly saving the purification run time and increasing the purification efficiency while still maintaining good separation with minimal resolution change.

Figure 4 shows one example of Sepax SRT-10/10C columns in 21.2 x 400 mm size that can run under very fast

flow rate at 7 mL/min and 10 mL/min which shortens the purification time to below 20 minutes while still achieving high resolution with minimal change.


Figure 4. SRT-10 SEC-300 21.2 x 400 mm at fast flow rates



Columns: SRT-10 SEC-300 (10 μ m, 300 Å 21.2 x 400 mm)
 Mobile phase: 150 mM phosphate buffer, pH 7.0
 Flow rate: 7.0 mL/min, 10.0 mL/min
 Detection: UV 280 nm
 Temperature: Ambient (23 °C)
 Injection volume: 3 mL
 Sample: 10 mg/mL BSA 66 kD; 0.17 mg/mL Uracil 120 D

With 0.75, 1 and 1.5 mL/min flow rate, Sepax SRT-10/10C 10 x 300 mm size columns show minimal resolution change.

Figure 5. SRT-10 SEC-300 10 x 300 mm at different flow rates



Column: SRT-10 SEC-300 (10 μ m, 300Å, 10 x 300 mm),
 Flow rate: 0.75, 1, 1.5 mL/min,
 Sample: 500 μ L mixture of Thyroglobulin 5.1 mg/mL, BSA 5.3 mg/mL, Ribonuclease A 5.2 mg/mL


| ID | 4.6x300 mm | 7.8x300 mm | 10x300mm | 21.2x300mm | 30x300mm | 50x300mm |
|--------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| Column Volume | 5 mL | 14 mL | 24 mL | 105 mL | 212 mL | 490 mL |
| V-injection | 5-100 μ L | 15-250 μ L | 25-940 μ L | 0.1-4.2 mL | 0.2-8.5 mL | 0.5-20 mL |
| Column Loading Guideline (BSA) | <= 1 mg | <= 2.5 mg | <= 9 mg | <= 42 mg | <= 85 mg | <= 235 mg |
| Standard Flow rate (Maximum) | 0.35 mL/min | 1.0 mL/min | 1.65 mL/min (2.0 mL/min) | 7.5 mL/min (10mL/min) | 15 mL/min (25 mL/min) | 41 mL/min (60 mL/min) |
| Typical Run Time | < 20 min | < 20 min | < 20 min | < 20 min | < 20 min | < 20 min |
| Back pressure | < 400 psi < 30 bar < 3 MPa | < 700 psi < 50 bar < 5 MPa | < 900 psi < 60 bar < 6 MPa | < 900 psi < 60 bar < 6 MPa | < 900 psi < 60 bar < 6 MPa | < 900 psi < 60 bar < 6 MPa |

Resolution

Benefiting from the rigid silica particle with high pore volume and uniform hydrophilic coating, SRT-10/10C packings provide a seamless transition from analytical separation scale up to preparative purification with high efficiency and high resolution for biomolecules.


Particle Size Impact on Resolution

Figure 6. The particle size impact on separation efficiency



Columns: SRT-10 10 μm, SRT 5 μm and Zenix 3 μm, SEC-300 (300 Å, 7.8 x 300 mm)
 Mobile phase: 150 mM phosphate buffer, pH 7.0
 Flow rate: 1.0 mL/min
 Detection: UV 214 nm
 Temperature: Ambient (23 °C)
 Injection volume: 5.0 μL
 Sample: 1) Thyroglobulin (670 kD); 2) BSA dimer (132 kD); 3) BSA (66 kD); 4) Ribonuclease A (13.7 kD); 5) Uracil (120 D).

Figure 7. The particle size impact on resolution - 21.2 x 300 mm




Columns: SRT-10 10 μm, SRT 5 μm, SEC-300 (300 Å, 21.2 x 300 mm)
 Mobile phase: 150 mM phosphate buffer, pH 7.0
 Flow rate: 7.0 mL/min
 Detection: UV 214 nm
 Temperature: Ambient (23 °C)
 Injection volume: 5.0 μL
 Sample: 1) Thyroglobulin (670 kD); 2) BSA dimer (132 kD); 3) BSA (66 kD); 4) Ribonuclease A (13.7 kD); 5) Uracil (120 D).

Column Length Impact on Resolution

The separation resolution can be greatly enhanced by increasing the column length. Figure 8 shows the protein separation by a 60 cm long SRT-10 SEC-300 column vs a 30 cm SRT-10 SEC-300 column. The efficiency with the 60 cm column is almost doubled that of the 30 cm column.

Figure 8. The column length impact on separation efficiency



Columns: SRT-10 SEC-300 (10 μm, 300 Å)
 Mobile phase: 150 mM phosphate buffer, pH 7.0
 Flow rate: 1.0 mL/min
 Detection: UV 214 nm
 Temperature: Ambient (23 °C)
 Injection volume: 5.0 μL
 Sample: 1) Thyroglobulin (670 kD); 2) BSA dimer (132 kD); 3) BSA (66 kD); 4) Ribonuclease A (13.7 kD); 5) Uracil (120 D).

Plate counts for SRT-10 SEC-300 columns with different lengths


| Peak | Proteins | 7.8 x 300 mm | 7.8 x 600 mm |
|------|----------------|--------------|--------------|
| 1 | Thyroglobulin | 691 | 1482 |
| 2 | BSA dimer | 2318 | 3452 |
| 3 | BSA | 3040 | 5827 |
| 4 | Ribonuclease A | 5272 | 10151 |
| 5 | Uracil | 14632 | 26821 |

Column ID Impact on Separation Resolution

Well packed SRT-10/10C semi-preparative and preparative columns also increase the separation efficiency and resolution in comparison to the analytical column with the same length.

Figure 9 shows the direct comparison of a 21.2 mm ID semi-prep column and a 7.8 mm ID analytical column, indicating that the semi-prep column increases the efficiency by 15%. The plate number of BSA is 4930 and 4338 for 21.2 x 300 mm and 7.8 x 300 mm columns respectively. For ribonuclease A, the plate number is 8098 and 6984 for 21.2 x 300 mm and 7.8 x 300 mm columns respectively.

Figure 9. The column ID impact on separation efficiency




Columns: SRT-10 SEC-500 (10 μ m, 500 \AA)
 Mobile phase: 150 mM phosphate buffer, pH 7.0
 Flow rate: 1.0 mL/min for 7.8 x 300 mm
 7.0 mL/min for 21.2 x 300 mm
 Detection: UV 214 nm
 Temperature: Ambient (23 $^{\circ}$ C)
 Injection volume: 5.0 and 20 μ L for 7.8 and 21.2 mm ID
 Sample: 1) Thyroglobulin (1.0 mg/mL), 670 kD; 2) BSA dimer, 132 kD; 3) BSA (1.0 mg/mL), 66 kD; 4) Ribonuclease A (1.0 mg/mL), 13.7 kD; 5) Uracil (2.5 μ g/mL), 120 D.

High Loading Capacity

Loading capacity is critical for size exclusion separation and purification. The following figures show high protein loading capacity. For example, more than 1 mg BSA can be loaded onto a SRT-10 SEC-300 7.8 x 300 mm analytical column. (The instrument for the following applications is AKTA Explorer FPLC)


Column Size: 7.8 x 300 mm

Figure 10. BSA (10 mg/mL) loading on SRT-10 SEC-300



Column: SRT-10 SEC-300 (10 μ m, 300 \AA , 7.8 x 300 mm)
 Mobile phase: 150 mM phosphate buffer, pH 7.0
 Flow rate: 1 mL/min
 Detection: UV 280 nm
 Temperature: Ambient (23 $^{\circ}$ C)
 Sample: 10 mg/mL BSA


Figure 11. QC protein standards loading on SRT-10 SEC-300



Column: SRT-10 SEC-300 (10 μ m, 300 \AA , 7.8 x 300 mm)
 Mobile phase: 150 mM phosphate buffer, pH 7.0
 Flow rate: 1 mL/min
 Detection: UV 280 nm
 Temperature: Ambient (23 $^{\circ}$ C)
 Sample: 5.5 mg/mL Thyroglobulin, 6.3 mg/mL BSA, 5.9 mg/mL Ribonuclease A


Column Size: 10 x 300 mm

Figure 12. BSA (10 mg/mL) loading on SRT-10 SEC-300



Column: SRT-10 SEC-300 (10 μ m, 300 \AA , 10 x 300 mm)
 Mobile phase: 150 mM phosphate buffer, pH 7.0
 Flow rate: 1-1.5 mL/min
 Detection: UV 280 nm
 Temperature: Ambient (23 $^{\circ}$ C)
 Sample: 10 mg/mL BSA


Figure 13. QC protein standards loading on SRT-10 SEC-300



Column: SRT-10 SEC-300 (10 μ m, 300 \AA , 10 x 300 mm)
 Mobile phase: 150 mM phosphate buffer, pH 7.0
 Flow rate: 1-1.5 mL/min
 Detection: UV 280 nm
 Temperature: Ambient (23 $^{\circ}$ C)
 Sample: Thyroglobulin 5.1 mg/mL, BSA 5.3 mg/mL, Ribonuclease A 5.2 mg/mL


Column Size: 21.2 x 300 mm and 400mm

Figure 14. High loading with 10 mg/mL BSA separation
SRT-10 SEC-300




Column: SRT-10 SEC-300 (10 μ m, 300 \AA , 21.2 x 300 mm)
 Mobile phase: 150 mM PB, pH 7.0
 Flow rate: 7 mL/min
 Detection: UV 280 nm
 Temperature: Ambient (23 $^{\circ}$ C)
 Sample: 10 mg/mL BSA

Figure 15. QC protein standards loading on SRT-10 SEC-300



Column: SRT-10 SEC-300 (10 μ m, 300 \AA , 21.2 x 300 mm)
 Mobile phase: 150 mM phosphate buffer, pH 7.0
 Flow rate: 7 mL/min
 Detection: UV 280 nm
 Temperature: Ambient (23 $^{\circ}$ C)
 Sample: 5 mg/mL Thyroglobulin, 5.4 mg/mL BSA, 5.22 mg/mL Ribonuclease A

Figure 16. SEC prep column- BSA loading study



Column: SRT-10 SEC-300 (10 μ m, 300 \AA , 21.2 x 400 mm)
 Mobile phase: 150 mM phosphate buffer, pH 7.0
 Flow rate: 7 mL/min, 17 bar
 Detection: UV 280 nm
 Temperature: Ambient (23 $^{\circ}$ C)
 Samples: 3 mL BSA (10, 20, or 30 mg/mL)

SEC prep column-BSA loading (See above chromatogram)


| BSA30mg | | | | | | | |
|----------------|---------|----------------|--------------|-------------|------|--------------------|-----------|
| Peak | RT(min) | Area (mAU*min) | Height (mAU) | W 1/2 (min) | Rs | Plates/meter (N/m) | Asymmetry |
| BSA dimer | 11.85 | 101.6082 | 102.906 | 0.8 | | 3069 | 0.34 |
| BSA | 13.07 | 417.1044 | 690.324 | 0.55 | 1.06 | 7720 | 0.88 |
| BSA60mg | | | | | | | |
| Peak | RT(min) | Area (mAU*min) | Height (mAU) | W 1/2 (min) | Rs | Plates/meter (N/m) | Asymmetry |
| BSA dimer | 11.76 | 195.445 | 190.721 | 0.87 | | 2535 | 0.34 |
| BSA | 13.02 | 1163.3833 | 3033.566 | 0.26 | 1.31 | 35524 | 0.86 |
| BSA90mg | | | | | | | |
| Peak | RT(min) | Area (mAU*min) | Height (mAU) | W 1/2 (min) | Rs | Plates/meter (N/m) | Asymmetry |
| BSA dimer | 11.82 | 28.065 | 20.617 | 0.86 | | 2640 | 0.33 |
| BSA | 13.19 | 2629.0897 | 4781.054 | 0.47 | 1.22 | 11005 | 0.64 |

Column: SRT-10 SEC-300 (10 μ m, 300 \AA , 21.2 x 400 mm)
 Mobile phase: 150 mM phosphate buffer, pH 7.0
 Flow rate: 7 mL/min, 17 bar
 Detection: UV 280 nm
 Temperature: Ambient (23 oC)
 Samples: 3 mL BSA (10, 20, or 30 mg/mL)

Comparison

The Sepax SRT-10 SEC column has better separation between Thyroglobulin and BSA (circled area). If lower molecular weight protein separation (between BSA and Ribonuclease A region, peak 3 and peak 4) is of interest, a smaller pore size SEC column is recommended. Sepax columns perform better with higher loading, better resolution, and shorter run time.

Figure 17. GE and Sepax 10 x 300 mm comparison



Column: SRT-10 SEC-300 (10 μ m, 300 \AA , 10 x 300 mm)
GE Superdex 200 (10 x 300 mm)

Mobile phase: 150 mM phosphate buffer, pH 7.0

Flow rate: 0.75 mL/min


Detection: UV 280 nm

Temperature: Ambient (23 $^{\circ}$ C)

Injection volume: 500 μ L

Samples: 5 mg/mL Thyroglobulin, BSA, Ribonuclease A

Figure 18. SRT-10 SEC-300 and Vendor S SEC comparison



Column A: SRT-10 SEC-300 (10 μ m, 300 \AA , 21.2 x 300 mm)

Flow rate: 2 mL/min

Injection volume: 200 μ L

Column B: Vendor S SEC (5 μ m, 20 x 300 mm)

Flow rate: 1 mL/min,

Injection volume: 100 μ L

Mobile phase: 50 mM Tris-HCl, pH 7.5, 150 mM NaCl, 10% glycerol


Detection: UV280 nm

Temperature: Ambient (23 $^{\circ}$ C)

Monoclonal Antibody Separation


Figure 19 shows that separation resolution of monoclonal antibody remains constant with increasing mAb loading.

Figure 19. Monoclonal antibody separation at various loadings




Sepax's SRT-10 SEC-300 gives consistently good resolution at different flow rates and is able to be run at higher flow rates than the GE Superdex 200.

Figure 20. Monoclonal antibody separation at different flow rates



Sepax's SRT-10 SEC-300 gives much sharper and taller peaks than GE's Superdex 200, which displays much broader peaks.

Figure 21. Comparison of SRT-10 and GE Superdex200



Column: SRT-10 SEC-300 (10 μ m, 300 \AA , 10 x 300 mm)
GE Superdex 200 (10 x 300 mm)

Mobile Phase: 150 mM Na₃PO₄, pH7

Flow rate: 0.75 mL/min

Instrument: AKTA Explorer FPLC

Detection: UV 280nm

Temperature: Ambient (23 $^{\circ}$ C)

Sample: 5 mg/mL MAb321

Injection volume: 200 μ L (1 mg)

High Robustness

SRT-10/10C packings have specially designed stationary phases that are densely bonded on the silica surface which enhances the stability of the column, resulting in high robustness at high flow rates.

High Stability

The proprietary stationary phases of SRT-10/10C packings utilize densely bonded chemistry on the silica surface, which greatly hinders the diffusion of the molecules, thus enabling high stability over a wide range of pH from 2 to 8.5. If a higher pH is being used, such as 9.0, equilibrate the column with 150 mM sodium phosphate at pH 7.0 and store the column in 150 mM sodium phosphate at pH 7.0.

Mobile Phase Compatibility

SRT-10/10C phases are compatible with most aqueous buffers, such as ammonium acetate, phosphate, trizma, etc. SRT-10/10C phases can tolerate high concentration of

salts, such as 2.0 M. Furthermore, SRT-10/10C columns are stable in both organic solvents (such as methanol, ethanol, THF, DMF and DMSO), and the mixture of water and organic solvents.

High Protein Recovery

SRT-10/10C phases are hydrophilic and neutral. Proteins and other biological molecules have negligible nonspecific interactions with SRT-10/10C stationary phases. The protein adsorption to the silica surface is suppressed, leading to a high recovery of intact proteins, which maintains the protein activity after separation. More than 95% recovery is achieved for BSA and lysozyme, the representatives for acidic and basic proteins, respectively.

Column Dimension Availability

SRT-10/10C SEC columns are available in dimensions 4.6, 7.8, 10, 21.2, 30 and 50 mm I.D., and 50, 100, 150, 250, 300 and 600 mm length. Sepax also offers custom-size columns.

SRT-10 Technical Specifications

| Phase | SRT-10/10C SEC-100 | SRT-10/10C SEC-300 | SRT-10/10C SEC-500 | SRT-10/10C SEC-1000 |
|--|--|--------------------|--------------------|---------------------|
| Material | Neutral, hydrophilic film bonded silica | | | |
| Particle size | 10 μm | 10 μm | 10 μm | 10 μm |
| Pore size (\AA) | ~100 | ~ 300 | ~ 500 | ~ 1000 |
| Protein MW range (native) | 100 - 30,000 | 5,000 - 1,250,000 | 15,000 - 5,000,000 | 50,000-7,500,000 |
| pH stability | 2 - 8.5 (pH 8.5-9.5 can be tolerated temporarily.) | | | |
| Salt concentration range | 20 mM - 2.0 M | 20 mM - 2.0 M | 20 mM - 2.0 M | 20 mM - 2.0 M |
| Maximum temperature ($^{\circ}\text{C}$) | ~ 80 | ~ 80 | ~ 80 | ~ 80 |
| Mobile phase compatibility | Aqueous and organic | | | |



Ordering Information

SRT-10 SEC Column

| ID x Length (mm) | SRT-10 SEC-100 (10 μm , 100 \AA) | SRT-10 SEC-300 (10 μm , 300 \AA) | SRT-10 SEC-500 (10 μm , 500 \AA) | SRT-10 SEC-1000 (10 μm , 1000 \AA) |
|------------------|--|--|--|--|
| | P/N | P/N | P/N | P/N |
| 50x250 | 225100-50025 | 225300-50025 | 225500-50025 | 225950-50025 |
| 30x300 | 225100-30030 | 225300-30030 | 225500-30030 | 225950-30030 |
| 30x250 | 225100-30025 | 225300-30025 | 225500-30025 | 225950-30025 |
| 21.2x300 | 225100-21230 | 225300-21230 | 225500-21230 | 225950-21230 |
| 21.1x250 | 225100-21225 | 225300-21225 | 225500-21225 | 225950-21225 |
| 21.2x50 | 225100-21205 | 225300-21205 | 225500-21205 | 225950-21205 |
| 10x300 | 225100-10030 | 225300-10030 | 225500-10030 | 225950-10030 |
| 10x250 | 225100-10025 | 225300-10025 | 225500-10025 | 225950-10025 |
| 10x50 | 225100-10005 | 225300-10005 | 225500-10005 | 225950-10005 |
| 7.8x300 | 225100-7830 | 225300-7830 | 225500-7830 | 225950-7830 |
| 4.6x300 | 225100-4630 | 225300-4630 | 225500-4630 | 225950-4630 |

SRT-10C SEC Column

| ID x Length (mm) | SRT-10C SEC-100 (10 μm , 100 \AA) | SRT-10C SEC-300 (10 μm , 300 \AA) | SRT-10C SEC-500 (10 μm , 500 \AA) | SRT-10C SEC-1000 (10 μm , 1000 \AA) |
|------------------|---|---|---|---|
| | P/N | P/N | P/N | P/N |
| 50x250 | 239100-50025 | 239300-50025 | 239500-50025 | 239950-50025 |
| 30x300 | 239100-30030 | 239300-30030 | 239500-30030 | 239950-30030 |
| 30x250 | 239100-30025 | 239300-30025 | 239500-30025 | 239950-30025 |
| 21.2x300 | 239100-21230 | 239300-21230 | 239500-21230 | 239950-21230 |
| 21.2x250 | 239100-21225 | 239300-21225 | 239500-21225 | 239950-21225 |
| 21.2x50 | 239100-21205 | 239300-21205 | 239500-21205 | 239950-21205 |
| 10x300 | 239100-10030 | 239300-10030 | 239500-10030 | 239950-10030 |
| 10x250 | 239100-10025 | 239300-10025 | 239500-10025 | 239950-10025 |
| 10x50 | 239100-10005 | 239300-10005 | 239500-10005 | 239950-10005 |
| 7.8x300 | 239100-7830 | 239300-7830 | 239500-7830 | 239950-7830 |
| 4.6x300 | 239100-4630 | 239300-4630 | 239500-4630 | 239950-4630 |

SRT-10/10C SEC Bulk Media

| SEC Phase | (10 μm , 100 \AA) | (10 μm , 300 \AA) | (10 μm , 500 \AA) | (10 μm , 1000 \AA) |
|-------------|--|--|--|---|
| | P/N | P/N | P/N | P/N |
| SRT-10 SEC | 225100-0000 | 225300-0000 | 225500-0000 | 225950-0000 |
| SRT-10C SEC | 239100-0000 | 239300-0000 | 239500-0000 | 239950-0000 |

How to Order

Please contact Sepax Sales Department:

Phone: (302)366-1101

Toll Free: 1-877-SEPAX-US

Fax: (302)366-1151

Email: sales@sepax-tech.com

5 Innovation Way, Suite 100
Delaware Technology Park
Newark, Delaware, 19711, USA

Discounts

Sepax Technologies offers the best discounts determined by the volume of the purchase. Please contact the Sepax Sales Department for your maximum discount.

Opening a Sepax Account

Call the Sepax Sales Department and supply your business information, and billing and shipping address to set up a Sepax account. Open account terms are subject to credit approval.

Payment Term

Terms of payment are net 30 days. Mastercard®, Visa®, and American Express® are accepted. There is no minimum order.

Return Policy

Shipping

If items are damaged in transit, simply follow these instructions:

- If shipment is visibly damaged on arrival, do not accept it until the delivery person has endorsed it with a statement for the extent of the damage.
- Notify us immediately of the damaged shipment in order for us to make the appropriate adjustment and/or provide you with return instructions.

Returns

- Sepax accepts eligible returns within 15 days of the customer receiving order.
- Non-eligible returns include products contaminated, treated, or tested, with isotope, radioactive chemical, or any other types of hazardous material, semi-prep and prep columns, custom products, bulk resins/materials, and demo purchase.
- Prior authorization is required for all returns. Please contact your local sales manager for prior authorization and Return Authorization Number.
- A 15% restocking charge will be made on all returns.
- Shipping costs are non-refundable. The customer pays for all of the shipping related costs including sending return product back to Sepax. A refund will only be processed upon receipt of the returned product.
- The return and refund will be made using the same method of purchase, i.e. through a distributor if purchased through distributor.

Warranty

Sepax Technologies warrants its products to be free from manufacturing defects for 90 days after the shipment. Sepax will accept for return or replacement any product which fails to meet the stated specifications. This warranty shall not apply to any defect, failure or damage caused by improper use or improper or inadequate maintenance and care. This warranty is exclusive and no other warranty, whether written or oral is expressed or implied. Sepax specifically disclaims the implied warranties of merchantability and fitness for a particular purpose. Under no circumstances shall Sepax be liable for direct, indirect or consequential damages arising from the use of its products. The maximum liability that Sepax will assume should be no more than the invoice price of the product.

*Please visit our website for the most updated literatures

Better Surface Chemistry for Better Separation



Sepax Technologies, Inc.

Delaware Technology Park
5-100 Innovation Way,
Newark, Delaware 19711

Tel: 302-366-1101
Fax: 302-366-1151

www.sepax-tech.com

