

# EZcount™ XTT Cell Assay Kit

**Product Code: CCK015**

## 1. Introduction:

Cell proliferation and death are essential processes for tissue generation and regeneration, organ development etc. in mammals and are usually under stringent control of extra and intracellular factors. Non-physiological alterations in levels of these factors lead to anomalous cytogenetic behavior of cells which in turn leads to cell transformation, uncontrolled cell growth - the initiating event for cancer development. Pharmaceutical research is hence largely focused on effects of drugs, cytotoxic agents and biologically active compounds which affect cytogenetics.

Multiple procedures are available for determination of cell proliferation and cytotoxicity. Simple and cheap methods for estimating cell viability (or death) are Trypan Blue exclusion and Erythrocin B staining. However, these methods are not sensitive enough and cannot be used for high throughput screening. Measuring the uptake of radioactive substances, usually tritium-labeled thymidine, is accurate but it is also time-consuming and involves handling of radioactive substances. Tetrazolium salts have been used to develop a quantitative colorimetric assay to estimate mammalian cell survival and proliferation. The assay detects living, but not dead cells and the signal generated is dependent on the metabolic state of the cells. This method can therefore be used to measure cytotoxicity, proliferation or activation. The results can be read on a multi-well scanning spectrophotometer or a standard ELISA reader and show a high degree of precision.

## 2. About the Assay:

The EZcount™ XTT Cell Assay kit is designed for determination of cell viability and effect of cytotoxic agent. This kit is based on quantitative measurement of extracellular reduction of XTT to water soluble orange coloured formazan derivative by metabolically active

cells. This reduction is mediated by mitochondrial dehydrogenase enzyme. Reduction of XTT is greatly enhanced by electron coupling reagent. Electron coupling reagent mediates reduction of XTT by transferring electron from cell surface/ cell membrane to extracellular XTT.

## 3. Applications:

- **Cell proliferation:** Quantification of changes in proliferative activity of cells caused by trophic factors, cytokines, and growth promoters.
- **Cell cytotoxicity:** Evaluation of effects of inhibitors or inducers of apoptosis, cytotoxic reagents, carcinogens and toxins.
- **Drug discovery:** High-throughput screening of various anti-cancer drugs.

## 4. Kit contents:

Contents		Kit Code	Storage
Code	Description	CCK003-1000*	
CCK015 Part A	XTT reagent	10 x 5ml*	-20°C in dark
CCK015 Part B	Electron coupling reagent	2 x 0.5ml*	-20°C in dark

\*Quantities supplied in excess to compensate operational losses

## 5. Materials required but not provided in the kit:

- Cells in appropriate medium without phenol red
- Adjustable pipettes and a repeat pipettor
- 96-well plate for culturing the cells
- 96-well plate reader capable of measuring the absorbance

## 6. General guidelines:

It is important to optimize experimental factors like cell density, incubation time, media composition and concentration of the agents under investigation prior to use of EZcount™ XTT Cell Assay Kit.

### Assay controls

- Include appropriate assay controls i.e.
  1. Medium control (medium without cells)
  2. Cell control (medium with cells but without the experimental drug/ compound)
  3. Vehicle control (medium containing the experimental drug or compound but no cells)
- Extracellular reducing components such as ascorbic acid, cholesterol, alpha-tocopherol, dithiothreitol present in the culture media may reduce the XTT to orange coloured formazan derivative. To account for this reduction, it is important to use the same medium in control as well as test wells.

### Accuracy:

- To obtain statistically significant data, perform the assay in triplicates or more.
- Accuracy of the assay depends on pipetting skills of the personnel. Inappropriate addition and mixing practices may result in erroneous and false-positive or false-negative results.
- Use of a repeating pipettor is recommended to deliver the reagents to the wells. This saves time and helps maintain more precise incubation times.
- Pipette tip should be equilibrated with the reagent before use. This is carried out by slowly filling up the tip with reagent and gently expelling the contents, several times.

### Incubation period:

- Different cell lines may have different properties such as metabolic activity and doubling time and hence respond to XTT differently. For this reason plating density and incubation period for every cell line should be optimized to obtain the results in linear range.

### Culture Medium

- Phenol red interferes with the measurement of orange colour; therefore the cell culture media used for this assay should not contain phenol red.

### Temperature

- Temperature affects the performance of the assay because of its effect on enzymatic rates. It is crucial to run the assay at a uniform temperature to ensure reproducibility across a single plate or among stacks

of several plates. Since absorbance or fluorescence reading are measured at room temperature, it is important to ensure adequate equilibration of assay plates after removal from a 37°C incubator to avoid differential temperature gradients. Stacking large numbers of assay plates in close proximity should be avoided to ensure complete temperature equilibration.

### Measurement of absorbance

- Absorbance can be read with a filter in the wavelength range of 450-500nm (primary wavelength).
- Reference wavelength (differential wavelength) should be between 630-690nm.

## 7. Directions for use:

***Users are advised to review entire procedure before starting the assay***

### 7.1 Preparation of XTT reagent:

1. Thaw XTT reagent and electron coupling reagent prior to use.
2. Mix 0.1ml of electron coupling reagent with 5ml of XTT properly to form activated XTT solution.
3. This amount will be sufficient for one 96 well plate (50µl/well).
4. Equilibrate the activated XTT reagent at room temperature before use.

### 7.2 Preparation of cells:

Always use freshly harvested cells for assay. Seed the cells in a cell culture flask or dish in an amount appropriate for the assay and incubate at 37°C in a 5% CO<sub>2</sub> environment. Allow the cells to grow up to 24 hours or till confluence is reached. Harvest the cells and use for the assay.

*(Note: Quantity of the cell suspension to be seeded in the medium depends upon doubling time of individual cell lines and seeding density to be used in assay).*

### 7.3 Procedure for determination of cell number or optimum cell density to be used in the assay:

Use the procedure given below to determine plating density and incubation period for cell line.

1. Harvest the cells as explained in section 7.2.
2. Adjust the cell density to  $1 \times 10^6$  cells/ml.
3. Serially dilute the cell suspension from  $1 \times 10^6$  to  $1 \times 10^3$  cells/ml using appropriate culture medium.
4. Seed 100µl of each dilution in 96-well plate in triplicate.

5. Include medium control in triplicate.
6. Incubate the cells under appropriate conditions depending on specific requirements for that cell line.
7. Add 50µl of activated XTT reagent to each well, including control.
8. Return the plate to incubator for 2 to 4 hours.
9. Once the incubation is over stir gently on a gyratory shaker to enhance uniform dispersion of reagent.
10. Read the absorbance on spectrophotometer or an ELISA reader by using 450nm as primary filter and 630nm as differential filter.
11. Determine the average values from triplicate readings and subtract from this value the average value for blank (i.e. medium control). Plot absorbance against cell density.
12. Number of cells to be used in the cell proliferation assay should lie within linear portion of the plot.

## 7.4 Assay procedures

### 7.4.1 Procedure for determining cell proliferation

1. Seed 100µl cell suspension in a 96-well plate at the required cell density, with or without the cell growth modifying agent.

*(Note:*

*a) If the cell growth modifying agent is a cytokine, metabolite, growth factor or any other compound, add its required quantity in the culture system.*

*b) If the cell growth modifying agent is any kind of radiation or waves, treat the cells with them for required period of time.)*

2. Incubate the plate at 37°C in a 5% CO<sub>2</sub> atmosphere for required period of time.
3. After the incubation period, remove the plates from incubator and add 50µl activated XTT reagent.
4. Return the plates to the incubator and incubate for 2 to 4 hours.

*(Note: Incubation time varies for different cell lines. Incubation time should be kept constant while making comparisons.)*

5. Remove the plate from incubator after incubation.
6. Stir gently on the gyratory shaker to enhance uniform dispersion of reagent.
7. Read the absorbance on spectrophotometer or an ELISA reader by using 450nm as primary filter and 630nm as differential filter.

### 7.4.2 Procedure for determining cell cytotoxicity

1. a) For adherent cells:  
Seed 100µl cell suspension in a 96-well plate at the required cell density, with or without the cell growth modifying agent. Allow the cells to adhere to the culture plate for about 24 hours. Add appropriate concentrations of the test reagent.
- b) For suspension cells:  
Seed 100µl cell suspension in a 96-well plate at the required cell density. Add appropriate concentrations of the test reagent immediately.
2. Incubate the plate for required period at 37°C in a 5% CO<sub>2</sub> atmosphere.
3. After the incubation period, remove the plates from incubator and add 50µl activated XTT reagent.
4. Return the plates to the incubator and incubate for 2 to 4 hours.  
*(Note: Incubation time varies for different cell lines. Incubation time should be kept constant while making comparisons.)*
5. Remove the plate from incubator after incubation.
6. Stir gently on the gyratory shaker to enhance uniform dispersion of reagent.
7. Read the absorbance on spectrophotometer or an ELISA reader by using 450nm as primary filter and 630nm as differential filter.

## 8. Storage and shelf life:


- XTT reagent and electron coupling reagent are light sensitive. Store them at -20°C in dark. Shelf life is 9 months.
- Use before expiry date given on the label.  
*Note: Repeated freezing and thawing should be avoided. Once thawed, remaining solution should be aseptically dispensed into smaller aliquots for further use.*

## 9. Processing of the data:

- Subtract the average 450nm absorbance values of the control wells from the average 450nm absorbance values of corresponding experimental wells.
- Plot 450nm absorbance versus concentration of test agent.
- For e.g. If A<sub>1</sub>, A<sub>2</sub> and A<sub>3</sub> are absorbances of test wells in triplicates at 570nm, the average absorbance A is,  
 $A = (A_1 + A_2 + A_3)/3$
- Similarly, if C<sub>1</sub>, C<sub>2</sub> and C<sub>3</sub> are absorbances of control wells in triplicates at 570nm, the average absorbance C is,  $C = (C_1 + C_2 + C_3)/3$

- Calculate normalized absorbance by subtracting C from A. Plot the values of normalized absorbances against concentration of test agent.

## 10. Performance characteristics:



The sensitivity of XTT to detect changes in cell number has been determined by plotting the graph of normalized absorbance values versus cell number.

CHO cells were serially diluted and treated with XTT reagent provided in EZcount™ XTT Cell Assay Kit, in a 96-well plate. After incubation for 0, 2 and 4 hours in a humidified incubator at 37°C, 5% CO<sub>2</sub>, absorbance was read at 450nm using an ELISA plate reader. The absorbance data was processed as given in point (9).

## 11. Advantages:

- **Time saving:** Absorbance can be measured directly without involvement of solubilization step.
- **Easy reagent preparation:** Ready to mix and pre-weighed reagents offer ease of reagent preparation
- **Reproducibility:** Entire assay can be performed in a single plate. Cells and reagents need not be transferred. This facilitates reproducibility of the results.
- **Sensitive and accurate:** Tetrazolium salt reduction strongly correlates with the metabolic activity of the cells. This allows use of very low cell densities.
- **Safety:** No radioisotopes are involved.
- **Fast:** Use of multi-well ELISA plates allows the processing of large number of samples.
- **Flexibility:** XTT works on adherent as well as suspension cell lines.

## 12. Troubleshooting guide:

Use the following troubleshooting guidelines for technical assistance

Problem	Cause	Solution
Coloration in XTT reagent	Microbial contamination or contamination with reducing agent	Discard the contaminated vial of XTT and prepare fresh reagent aseptically
	Excessive exposure of reagent to light	Store in dark at -20°C
Very high absorbance values	Too much reduction of XTT due to high cell densities	Repeat the assay with reduced cell densities
	Too much reduction of XTT due to long incubation period	Repeat the assay with reduced incubation period
	Microbial contamination	Discard. Repeat the assay with new media and reagents
Very low absorbance values	Very low cell density	Repeat the assay with high cell densities
	Small incubation period	Repeat the assay with longer incubation period. Certain cell types require longer incubation period of up to 24 hours
	Improper selection of filter for reading the absorbance	Choose appropriate filters within the range of 450-500nm
Random absorbance values/ poor consistency of replicates	Inefficient pipetting techniques	Perform the assay using automated electronic pipettes for seeding the cell suspension and adding the reagents
	Test compound under study which is responsible for improper response of the cells to XTT	Refer to the pharmacological properties of the compound
Blank/ medium control (i.e. medium without cells) give high absorbance readings	Microbial contamination	Discard. Repeat the assay with new media and reagents.

Revision No.: 1//2012

### Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt. Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.