

BIOPROCESS PRODUCTS



ENTERPRISE INTRODUCTION

Founded in 2016, Wuhan CEKG Technology CO., LTD. is a developer, manufacturer and marketer of high-tech medical equipment and consumables. With its core skills in cell culture technology, the company is committed to providing a range of high-quality products, e.g. traditional bioreactor, rocker bioreactor, single-use bioreactor, fabric carrier, single-use bags, single-use magnetic mixer, fabric carrier bag etc., as well as offering cell culture technology services.

The company has its own cell culture laboratory, which can provide cell culture technology services and bioreactor maintenance and certification services.

The company has successively applied for 33 patents, won the honorary title of high-tech enterprise and Wuhan science and technology "Little Giant" enterprise, and passed the ISO9001 quality system certification.

The company has its own main brand CEKG®, which focuses on providing cell culture equipment and consumables for vaccines, cell and gene therapy, antibodies and other industries. After the unremitting efforts of the team, the company has become an advanced high-tech enterprise in the industry.



QUALIFICATION HONOR







PATENTS
















CONTENTS

ABOUT CEKG

SINGLE-USE BAG	01-10
SKC400 SINGLE-USE MAGNETIC MIXING BAG	11-12
SKC400 SINGLE-USE MIXING SYSTEM	13-14
SKC400 ROCKER MIXER	15-16
SKC400 STORAGE SYSTEM	17-19
FABRIC CARRIER	20-23
SKC600 ROCKER BIOREACTOR	24-27
SKC600 200L ROCKER BIOREACTOR	28-29
SKC900 SINGLE-USE STIRRED BIOREACTOR	30-33
PROCESS APPLICATION	34-38
SKC300 GLASS BIOREACTOR	39-42
SKC800 STAINLESS STEEL BIOREACTOR	43-46
SKC700 SINGLE WALL BIOREACTOR	47-49
SERVICES	





PRODUCT INTRODUCTION

Low-Density Polyethylene (outer layer)

Bonding layer

Ethylene, vinyl alcohol copolymer

Bonding layer

Ultra-low density polyethylene (Liquid contact layer)

Ultra-low density polyethylene (outer layer)

Ultra-low density polyethylene

Bonding layer

Ethylene, vinyl alcohol copolymer

Bonding layer

Ultra-low density polyethylene

Ultra-low density polyethylene (Liquid contact layer)




SINGLE-USE BAG PRODUCTS


FILM INTRODUCTION

All single-use bags of CEKG are made with Multi-layer PE co-extruded film

FEATURE

















COMPATIBILITY

Items	Reference standards
Hemolysis test	ISO 10993-4
Cytotoxicity	ISO 10993-5
Implantation test	ISO 10993-6
Irritation and Sensitization test	ISO 10993-10
Acute systemic toxicity test	ISO 10993-11
Bacterial endotoxins-LAL test	USP<85>
Biological reactivity testing,in vivo, class VI	USP<88>
Plastic container European Pharmacopoeia test CH. 3.1.5	USP<661>

CHARACTERISTICS


Name	Test Value (pre-irradiation/ post-irradiation)	Reference Standards
Haze	7/7%	ASTM D-1003
Transparency	97/97%	ASTM D-1003
Light transmittance	93/93%	ASTM D-1003
Tensile strength	14/13MPa	ASTM D-882
Elongation at break	280/300%	ASTM D-882
Elasticity coefficient	370/350%	ASTM D-882
The lowest tolerable temperature	-45/-45°C	ISO8570
Density	0.9g	ASTM D-792
Water vapor permeability	0.35/0.32 (g/ (m².day)) (23°C,100%RH)	ASTM F-1249
Oxygen permeability	<0.05/<0.05 (cm³/ (m²)) (23°C,0%RH)	ASTM D-3985
Carbon dioxide permeability	<0.2/<0.2 (cm ³ / (m ²)) (23°C,0%RH)	ASTM D-2476








SINGLE-USE 2D BAG



FEATURE







03 Application for storage and transfer



Customization for configuration

ADVANTAGE

PRODUCING IN IS07/IS05 71

LEAKAGE TEST



Production in IS07/IS05

100% leakage test

FREEZED

73

CUSTOMIZED




Freezed

Customized

APPLICATION












SINGLE-USE 3D BAG



FEATURE



Good stability



Good compatibility



Application for storage and transfer



Customization for configuration

ADVANTAGE

PRODUCING IN IS07/IS05



LEAKAGE TEST



Production in IS07/IS05

100% leakage test

FREEZED



CUSTOMIZED



Freezed

Customized

APPLICATION



Storage and transfer of medium



Storage and transfer of buffer



Collection and storage of samples and compositions




Storage of intermediate products and semi-products








SINGLE-USE MAGNETIC MIXING BAG



FEATURE







Customization for configuration

APPLICATION

(01) Medium preparation

Buffer preparation

pH adjustment

Intermediate products and semi-products preparation

ADVANTAGE

PRODUCING IN IS07/IS05

LEAKAGE TEST

Production in IS07/IS05

100% leakage test

FREEZED

CUSTOMIZED



Freezed


Customized

(05) Intermediate products mixing

Intermediate products inactivation

Semi-products preparation

SINGLE-USE SAMPLING BAG



FEATURE



Good chemical compatibility and biocompatibility



Customized to 5 sampling bag kits/10 sampling bag kits



3 ports and 2 ports optional, and easy for configuration to suit sampling process



Assembling with 2D bags/3D bags/mixing bags/cell culture bags

APPLICATION



Sampling in various Processes

ADVANTAGE

PRODUCING IN IS07/IS05



LEAKAGE TEST



Production in IS07/IS05

100% leakage test

FREEZED



CUSTOMIZED



Freezed

Customized





CONSUMABLES

CPC Connectors

Different connectors including quick connectors, aseptic connectors, aseptic dis-connectors and SIP connectors.



CPC quick connectors

Quick connectors including MPC.MPX, MPU and sanitary connectors; sanitary connectors can be connected to sanitary clamps with integrated into single-use or stainless steel process system





\

CPC aseptic connectors

CPC aseptic connectors can be applied into no-clean area and various types optional.





03

CPC aseptic dis-connectors

CPC aseptic dis-connectors can avoid contamination and various types optional





04

CPC SIP connectors

CPC SIP connectors are applied under SIP, and various types optional







SINGLE-USE ADDITION BOTTLES/3 PORTS BOTTLE CAP

Single-use addition bottle

Pre-sterilized single-use bottle kit for ready to use C-flex tubing can be connected to cell culture bags

3 ports bottle cap

3 ports bottle caps match with glass bottles



SILICONE TUBING/C-FLEX TUBING

SILICONE tubing/C-flex tubing from domestic or abroad brands for bio-pharmaceutical process





CONNECTORS AND TUBING CLAMPS

FEATURE



Including different sizes and types




Barbs specially designed for silicone, C-flex and other hoses, easy to install and uneasy to pull out




Durable and smooth inner surface



For piping connection/reducing conversion/ branches or collection in bio-pharmaceuticals









CUSTOMIZED SERVICE

SERVICE DESCRIPTION

Provide customization of single-use cell culture bags/assembling/irradiation services according to the demands of customersAccessories include hoses/quick connectors/aseptic connectors/SIP connectors/sampling bottles/filters etc.

FEATURE



Provide customization of single-use cell culture bags/process design



Production line process to ensure quality



Production in ISO5/ISO7




Fast and partial delivery

APPLICATION

01

Customization cell culture bags and assembly services

Provide customization of singleuse cell culture bags/process design



02

Customization filter and piping assembly

Composed of one or more filters and piping and the kit sterilized



03

Customization bottle and piping assembly


Composed of one or more bottles and piping and the kit sterilized





Customization bags and assembly

Composed of more sampling bags/ addition bags/batch bags and the kit sterilized





Customization bottles and assembly

Composed of addition bottles/ harvest bottles/batch bottles and the kit sterilized











SKC400 SINGLE-USE MAGNETICMIXINGSYSTEM

Container Container Coperation door Load cell

APPLICATION



Medium mixing and prepared



Mixing for vaccine adjuvant



Mixing pH or conductivity of intermediate product in purification process



Ultra filtration and semifinished product prepared



Dilution preparation of concentrated buffer

INTRODUCTION

SKC400 single-use magnetic mixing system adopts brand components from domestic and abroad. The equipment has a neat appearance, easy to clean and disinfect, concise human-machine interface, stable process control and convenient operation.

SKC400 single-use magnetic mixing system is a stable and high efficient system with stainless steel liquid container to meet requirements of liquid-liquid/solid-liquid mixing etc.

SKC400 single-use magnetic mixing system composed of master controller, universal magnetic drive unit, single-use bag with magnetic impeller, and a stainless steel container.

FEATURE



Piping skid and container manufactured by SS 304



The magnetic impeller in bag driven by coupling motor without seal



Siemens PLC control, color LCD touch screen, control stability and easy to operate



External needle printer can print the real data optional The single-use mixing bags are independently produced by ourselves with validation and stable long-term supply.



3 level security, audit trail, batch report, electronic signature, and GMP standards.



Modular design and configuration with temperature/load cell/pH/ DO/conductivity, etc.

SPECIFICATION

Blade quantity	4/3 blades
Materials	PP
Maximum speed	300rpm
Location	Center or eccentricity





SPECIFICATION

Control System	Siemens color touch screen, Ethernet interface, USB interface, Human-machine interface, and save data
Agitator	Magnetic drive.Max. speed: 300rpm
Impeller*	Single-use magnetic impeller, impeller type optional. Impeller material is PP
Agitating bag ※	Standard bags or customized according to container size and volume
Container Type *	Round and square
Load cell*	Optional
рН፠	Optional
Conductivity*	Optional
Temperature*	Optional
DO	Optional
Jacket design*	Optional. for low-temperature mixing
Software design *	Optional. 3 level security, audit trail, batch report, electronic signature, and GMP standards
Needle printer*	Optional
Operating temperature	No more than 45°C

Note: *Optional

SPECIFICATION

Volume	50L	100L	200L	500L	1000L	2500L
Design	Discharge from	Discharge from bottom, sampling ports/addition ports and probe ports in the side				
Material	SS304	SS304				
Surface polished	Ra≤0.6μm					
Trundle	Large castor and polyurethane material					
Installation	Side operation door					







SKC400 SINGLE-USE MIXING SYSTEM



INTRODUCTION

It adopts brand components from domestic and abroad. The equipment has a neat appearance, easy to clean and disinfect, simple man-machine interface, and stable process control.

It is mixing system with advantages of portability/high mixing efficiency/convenient use/small footprint and convenient transportation.

Applied for liquid-liquid/solid-liquid mixing and meets the demands for small/medium scale mixing system.

FEATURE



Piping skid and container manufactured by SS 304



Transfer base



The upper stirring matched with the turbine stirring Impeller, and PTFE optional



Volume 50L/100L/200L optional



Siemens PLC control, color LCD touch screen, stability and easy to operate



Standard bags and customization



Emergency switch and 3 level security



Container SS304 and HDPE optional



Platform scale optional

APPLICATION



Mixing of medium and prepared



Mixing of intermediate product in purification process



Mixing of Ultra filtration process

SPECIFICATION

Control system	Servo mixing system, Siemens color touch screen, Ethernet interface, USB interface
Agitator	Upper mixing, Max. speed: 300rpm
Impeller*	Turbine agitator, blade material is PTFE
Agitating bag *	Standard bags or customization
Container type *	Round
Load cell ※	Optional
Operating temperature	No more than 45°C

Note: *Optional

SPECIFICATION

Volume	50L	100L	200L	
Material	SS304			
Surface polish	Ra≤0.6μm			
Trundle	Large castor and polyurethane material			
Bag installation	Easy to take out and install			
Drain	Bottom or drain optional			

Note: **⊗Optional**







SKC400 ROCKER MIXER

PRODUCT INTRODUCTION



INTRODUCTION

The rocker mixer adopts non-intrusive wave mixing, which avoids the shearing force damage to the mixture by the end of the stirring blade. Compared with the stirring type mixer, it can provide a mild and low-shear mixing condition.

For cells/viruses/plasmids and decomposable medium, our company has developed a low temperature(5°C or below) rocker mixer and meet the requirements of low shear force with high mixing efficiency.

FEATURE



Low-shear mixing condition



Customized for trays and bags optional



No other stirring blades, saving transfer time and cost



Modular design and configuration with temperature/load cell/pH/ DO/ conductivity, etc.



Standardized operation process and easy to verify

APPLICATION



Plasmid mixing and plasmid lysis



Mix the virus liquid under low shear force and low temperature



Dilution preparation of concentrated buffer



Mixing pH or conductivity of intermediate product



Part of ADC reaction container



Preparation of semi-finished products



Mixing in a low temperature condition or easy decomposed medium



Low temperature mixer

SPECIFICATION

Tray size	10L	20L	50L	200L
Total Mixing bag volume	2L/5L/10L	10L/20L	20L/50L	100L/200L
Processing capacity	1L - 5L	2-10L	5-25L	20-100L
Control system	Siemens color touc	h screen, Ethernet inte	rface, USB interface	
Mixing system	Mixing time, rocker	speed and rocker angl	e (if any) automatical	ly stop after mixing.
нмі	Alarm record interfa	ce/trend graph interface	e/calibration interface	/time cascade interface, etc.
Rocker speed	1-40rpm			1-25rpm
Lid ※	Optional			- 1 %
Tilt angle※	Optional. Range: 1-1	2°		1-10°
Conductivity detection*	Optional			
pH detection※	Optional. Range: 5-9			
Temperature control*	Optional. PID autom	atic control, accuracy: ±	: 0.2°C; Temperature co	ntrol range: 5-40°C max. 60°C
Low temperature control ::	Optional. Low tempe	erature refrigerant to cor	ntrol the low temperatu	ire, PID control, 5 °C or below
Pumps ※	Optional. 1-3 fixed-s	peed/variable-speed pe	ristaltic pumps	
Load cell ※	Optional			
SCADA Software **	Optional. 3 level sec	urity, audit trail, batch re	port, electronic signat	ure, and GMP standards
Power Supply	AC220V,50/60Hz,1	0A		AC220V, 50/60Hz, 16A

Note: *Optional







SKC400 STORAGE SYSTEM

PRODUCT INTRODUCTION

Load cell

Container



Operating Port

INTRODUCTION

SKC400 single-use storage system adopts brand components from domestic and abroad. The equipment has a neat appearance, easy to clean and disinfect, concise human-machine interface, stable process control and convenient operation.

SKC400 single-use storage system is a stable and high efficient system with stainless steel liquid container to meet requirements of low temperature storage.

SKC400 single-use storage system composed of master controller, single-use storage bags, and stainless steel containers.

Single-use storage bags are independently produced by ourselves with validation and stable long-term supply.

FEATURE



Piping skid and container manufactured by SS 304



Siemens PLC control, color LCD touch screen, stability and easy to operate



Modular design and configuration with temperature/load cell/pH/DO/conductivity, etc.



External needle printer can print the real data optional



The single-use storage bags are independently produced by ourselves with validation and stable long-term supply



3 level security, audit trail, batch report, electronic signature, and GMP standards



Customized design, 2/3 layer stacking, operating door, etc.



Unnecessary of CIP and SIP

APPLICATION



Medium storage



Buffer storage



Low temperature storage



Intermediate storage



Continuous low-temperature harvesting and storage of medium or virus medium

SPECIFICATION

Control system *	Siemens color touch screen, Ethernet interface, USB interface			
Storage bag *	Standard bags or customization			
Container Type ※	Round and square			
Load cell*	Optional			
рНЖ	Optional			
Conductivity*	Optional			
Temperature*	Optional			
Dissolved oxygen *	Optional			
Jacket design※	Optional. For low-temperature mixing			
Software design*	Optional. 3 level security, audit trail, batch report, electronic signature, and GMP standards			
Needle printer*	Optional			
Operating temperature	No more than 45°C			

Note: *Optional

SPECIFICATION

Volume	50L	100L	200L	500L	1000L	2500L
Design	Discharge from	n bottom. Sam	pling ports/ad	dition ports ar	nd probe ports	in the side
Material	SS304	SS304				
Surface polished	Ra≤0.6μm					
Trundle	Large castor and polyurethane material					
Installation	Side operation	Side operation door				





CONTAINER

PP material, used for stacking storage and handling of a large amount of liquid, foldable placement, saving space, 50L, 200L, 500L, 1000L, etc. optional.









HDPE material, used for the storage and transportation. 50L, 200L, 500L, 1000L, etc. optional.

SPECIFICATION

Volume	50L	100L	200L	
Material	HDPE plastic	HDPE plastic		
Trundle*	Large castor and polyurethane material			
Installation*	Side operation door			

Round container

Transfer Base

Note: ***Optional**

SPECIFICATION

volume	50L	100L	200L	500L	1000L	2500L
Material	SS304					
Installation*	Side operation door					



Note: *Optional

PROJECT SITE








PRODUCT INTRODUCTION

FABRIC CARRIER

Three-layer diamond fabric carrier

Disk fabric carrier



Three-layer diamond fabric carrier and disk fabric carrier, the adherence surface area is 1200 cm2/g.







DESCRIPTION

The fabric carrier is a carrier applied for the cultivation of animal cell adherent cells, which provides a large adherent area for cell growth and realizes high-density cell culture.

The fabric carrier is composed of stable polyester fiber and polypropylene stent suitable for cell growth. It is made by welding, molding, cleaning, charge treatment and other processes. Animal-derived component-free(ADCF)

APPLICATION


It is suitable for mammalian/insect cell culture. It is mainly used for the production of secreted products, e.g. genetic recombination proteins/virus vaccines etc. Now, it can also be used to harvest cells, e.g. stem cells/immune cells etc.

Fabric carriers have been widely used in R&D, e.g. human rabies vaccines/IPV/EV71 vaccines /influenza vaccines/COVID-19/swine fever vaccines/avian flu vaccines/pseudo-rabies vaccine/porcine diarrhea vaccines/EPO/TPO protein drugs/Viral vectors (LV, AAV, ADV, RV, HSV, etc.)/cell therapy, etc. The cell lines that have been cultured include: VERO, CHO, BHK, MDCK, ST, MARC-145, ES, MSC, CEF, MRC-5, SF9, HEK-293, 123A, 127A, GAMMA and primary cells.






STRUCTURE




VIEW








Electron micro-graph of fabric carrier image









Electron micro-graph of cell image

INFRARED SPECTROGRAM

Soak in ultrapure water overnight, 121°C, 60min sterilization



Soak in 0.5mol NaoH for 24H





ADVANTAGE

CERTIFICATION

21

Meet requirements of declarations in domestic and abroad.

HIGH DENSITY

73

High area/volume ratio with high density.

PERFUSION

22

Separate the medium/cells/harvest products/ perfusion culture.

EASY TO CELL CULTURE




Similar to the static state of a spinner flask. Simple and easy to control.

CELL CULTURE CASE

Wuhan CEKG Technology Co., Ltd. data sheet of 2L volume fabric carrier bag Cell type: VERO; Inoculation density: 3.5*10⁵ cells/ml; Working volume: 1L; Growth days: 5 days.

VERO cell growth curve



ORDER INFORMATION

Model NO.	Description (Name, Specification, Size, Configuration)
SKC0101-0050	Fabric carrier/Disk-Fabric, double-layer, 50g/bag
SKC0101-0050G	Fabric carrier/Disk-Fabric, double-layer, Gamma irradiation, sterility, 50g/bag
SKC0101-0250	Fabric carrier/Disk-Fabric, double-layer, 250g/barrel
SKC0101-0250G	Fabric carrier/Disk-Fabric, double-layer, Gamma irradiation, sterility, 250g/barrel
SKC0101-1000	Fabric carrier/Disk-Fabric, double-layer, 1000g/bag
SKC0101-1000G	Fabric carrier/Disk-Fabric, double-layer, Gamma irradiation, sterility, 1000g/bag
SKC0102-0050	Fabric carrier/Diamond-Fabric, three-layer, 50g/bag
SKC0102-0050G	Fabric carrier/Diamond-Fabric, three-layer, Gamma irradiation, sterility, 50g/bag
SKC0102-0250	Fabric carrier/Diamond-Fabric, three-layer, 250g/bag
SKC0102-0250G	Fabric carrier/Diamond-Fabric, three-layer, Gamma irradiation, sterility, 250g/bag
SKC0102-1000	Fabric carrier/Diamond-Fabric, three-layer, 1000g/bag
SKC0102-1000G	Fabric carrier/Diamond-Fabric, three-layer, Gamma irradiation, sterility, 1000g/ba



CRITERION


Test item	Specification	Test standard		
Bioburden	<0.5 (CFU/g)	ISO 11737-1:2006		
Endotoxin	<0.05 (EU/ml)	Chinese Pharmacopoeia 2015.Vol IV:1143		
Cytotoxicity	Grade"0"	ISO 10993-5-2009		


















SKC600 BIOREACTOR

PRODUCT INTRODUCTION



INTRODUCTION

SKC600 bioreactor is developed for cell culture. It uses a rocker platform to place single-use suspension cell culture bags or single-use fabric carrier culture bags(produced by ourselves). Also It can control the rocker speed, tilt angle, gas flow, temperature, pH, DO, concentration of CO2 and O2 and weight.

SKC600 bioreactor supports the replacement of trays of different sizes. 10L trays (2L, 5L, 10L culture bags), 20L trays (20L culture bags) and 50L trays (50L culture bags) optional.

CONTROL MODE

Independent control and tray replacement

One controller can control one rocker and replace different trays easily.



Parallel control and multiple control

2-8 rockers can be controlled by a computer.








Testing site

НМІ



FEATURE

Tray replaced rapidly
10L trays, 20L trays, 50L trays replaced easily

Suspension culture and adherent culture
Suitable for different cell culture bags, e.g. suspension cell culture and adherent cell culture. Continuous

perfusion culture for adherent cell is achieved easily with our patented fabric carrier bag

Batch culture and perfusion culture of suspension cells Perfusion membrane cell culture bag is used for the perfusion culture and high-density culture of suspension cells and the N-1 scale-up

Adherent cell trypsinization, transfer and scale up Fabric carrier cell culture bags can achieve amount cell culture adherent, and then trypsinizate and prepare amount adherent seed cells into next process

SPECIFICATION

Tray size		10L		20L	50L		
Total volume (suspension culture)	2L	5L	10L	20L	50L		
Total volume (adherent culture)	2L	5L	10L	20L	50L		
Max. working volume	1L	2.5L	5L	10L	25L		
Control system	Siemens PLC, 15-inch color touch screen, Ethernet interface, USB interface						
Temperature control	Heating plate, PID automatic control, accuracy: ±0.2°C; Temperature control range: 5-40°C above						
Rockerspeed	1-40rpm						
Tilt angle	1-12°						
Gas ※	3 thermal mass flow controllers (Air/O2/CO2), 1 thermal mass flow controller with 3 gas mixer optional						
Exhaust *	Exhaust condenser, exhaust filter heater optional						
Exhaust analysis	02(0-50%), CO2 (0-15%)						
Pressure detection	Range:-100mbar100mbar; Stop air intake when it is exceeded.						
рН	Single-use pH probe, control range: 6-8, accuracy: ±0.05pH; PID control; cascaded by gas flow, peristaltic pump, etc.						
DO	Single-use DO probe, control range: 0-100%, accuracy: ±3%; PID control; cascaded by rocker speed, MFC, etc.						
Pumps※	4 fixed-speed peristaltic pumps, customized feeding/acid feeding /alkali feeding/ discharging etc. 25RPM or 100RPM variable speed optional						
Load cell	Load: 0-50kg, accuracy: \pm (0.050 + 1% loading sampling) kg; cascaded by pumps						
Utilities	Air/O2/CO2, pressure: 15-20PSIG						
Signal interface*	Optional						
Power supply	AC220V, 50/60HZ, 10A						

Note: *Optional





SKC600 200L BIOREACTOR



200L INTRODUCTION

SKC600 200L bioreactor is developed for cell culture. It uses a rocker platform to place single-use suspension cell culture bags or single-use fabric carrier culture bags(produced by ourselves). Also It can control the rocking speed, tilt angle, gas flow, temperature, pH, DO, concentration of CO2 and O2 and weight.

It is used for GMP production line, and also prepare suspended seed cells or adherent seed cells for large-scale stainless steel/single-use bioreactors more than 200L.

One 200L or two 100L cell culture bags optional.

200L FEATURE



For suspension cell culture and adherent cell culture



Cell culture bags are produced by ourselves with validation and stable long-term supply



Compatible with batch cell culture and perfusion cell culture



Monitored the whole production process and software to meets requirements of GMP



Trypsinization, transfer and scale up



Freeze-thawing and lysis



Process scaled-up

200L APPLICATION



Cultivation of mammalian cells/insect cells/ plant cells/bacteria/immune cells/stem cells.



R&D and production for virus vaccines



Seed preparation for suspension and adherent cells.



R&D and production for protein products



Preparation of virus seeds for suspension and adherent cells.



R&D and production for cell growth factors



Suspension culture and adherent culture of CAR-T, CAR-NK and other immune cells in cell therapy.



R&D and production for viral vectors



Suspension culture and adherent culture for IPS and MSC etc.



R&D and production for cell exosomes

200L SPECIFICATION

Control system	Siemens PLC, 15-inch color touch screen, Ethernet interface, USB interface
Temperature control	Heating plate, PID automatic control, accuracy: ±0.2°C; Temperature control range: 5-40°C above
Rocker speed	1-25rpm
Tilt angle	1-10°
Gas ※	3 thermal mass flow controllers (Air/O2/CO2), 1 thermal mass flow meter with 3 gas mixer optional
Exhaust **	Exhaust condenser, exhaust filter heater optional
Exhaust analysis	02 (0-50%), CO2 (0-15%)
Pressure detection	Range: -100mbar: —-100mbar; Stop air intake when it is exceeded.
pН	Single-use pH probe, control range: 6-8, accuracy: ±0.05pH; PID control; cascaded by MFC, peristaltic pump, etc
DO	Single-use DO probe, control range: 0-100%, accuracy: ±3%; PID control; cascaded by rocker speed, MFC, etc.
Pumps*	4 fixed-speed peristaltic pumps, 25RPM or 100RPM variable speed optional
Load cell	Load: 0-200kg, accuracy: \pm (0.050 + 1% loading sampling) kg; cascaded by pumps optional
Utilities	Air/O2/CO2, pressure: 15-20PSIG
Signal interface※	Optional
Power supply	AC220V, 50/60HZ, 10A

Note: ***Optional**






SKC600 MATCHING BAG

PRODUCT INTRODUCTION







SUSPENSION CELL CULTURE BAG AND PERFUSION MEMBRANE CELL CULTURE BAG

INTRODUCTION

Suspension cell culture bag can culture suspension cells.

Perfusion membrane cell culture bag can carry on metabolism of suspension cell culture to achieve ultra-high density cell culture and continuous perfusion cell culture.



Addition bag

Discharge bag

FEATURE

N-1 SCALE UP



Scale-up of N-1 can be achieved easily, saving 1-2 larger bioreactors

STERILE AND ENCLOSED



Good stability, sealed condition and less contamination

HIGH DENSITY PERFUSION



Perfusion membrane cell culture bag

Ultra-high-density cell culture and continuous perfusion cell culture

CUSTOMIZED



Customized

APPLICATION



Seed cell culture and N-1 scale up cell culture of protein drugs/antibodies



Higher density cell culture of mammalian cells/insect cells/plant cells/bacteria/ immune cells/stem cells



Preparation of virus seeds



Preparation of virus vaccines, e.g. swine fever/rabies/flu etc.



Suspension cell culture and enrichment of CAR-T, CAR-NK etc.




Suspension cell culture for IPS and MSC etc.



FABRIC CARRIER CELL CULTURE BAG

INTRODUCTION

Fabric carrier bag can achieve adherent cells culture by suspension mode. Fabric carrier is pre-filled into single-use bag and pre-sterilization. The adherent cells will adhere and grow in the fabric carrier and the fabric carrier with cells suspend in the medium. The cells or secreted products can be harvested through ourselves' patented technology.



FEATURE

MATRIX AND REPLACEMENT

71

Large adherence area, instead of cell factory, bottles to prepare the seed cell and seed virus

TRYPSINIZATION & TRANSFERANCE

73

To achieve adherent cell culture scale up by trypsinization

SAMPLING FEBRIC CARRIER

2

Sampling of fabric carrier and counting cell density

CONTINUOUS PERFUSION

94

Continuous perfusion cell culture

HIGH DENSITY

75

High-density cells and viruses

CERTIFICATION

706

Meet requirements of declarations in domestic and abroad

STERILITY AND STABILITY

77

Good stability and single-use enclosed condition

FREEZE THAWING AND SPLITTING

98

Freeze-thawing and lysis of bags

APPLICATION



Adherence cell culture of mammalian cells/ insect cells/plant cells/bacteria/immune cells/ stem cells



High-density cell culture of hybridoma cells/ stem cells etc.



Seed and virus preparation for adherent cells.



R&D for cell growth factors



Adherent cell culture of IPS and MSC etc.



R&D and production process for viral vectors, e.g. LV (Lentivirus)/ADV(Adenovirus)/AAV (Adeno-associated virus)/RV (Retrovirus)/HSV (Oncolytic virus) etc.



R&D and production for protein drug, e.g. EPO/TPO/TPA, etc.



R&D and production process for cell exosomes



R&D and production for human virus vaccines, e.g. COVID-19 vaccine/rabies vaccine/EV71 vaccine/IPV vaccine/influenza vaccine/vaccinia, etc.



R&D and production for animal virus vaccines, e.g. swine fever/PRRS/swine diarrhea/swine gastroenteritis/avian influenza/pseudorabies etc.

COMPARASION DATA BETWEEN FABRIC CARRIER AND ALTERNATIVE PRODUCTS

Total volume of carrier bag	2L	5L	10L	20L	50L	100L	200L
Working volume	1L	2.5L	5L	10L	25L	50L	100L
Fabric carrier quantity	30g	75g	150g	300g	750g	1500g	3000g
Adherence area	36000cm ²	90000cm²	180000 cm ²	360000cm ²	900000cm ²	1800000 cm ²	3600000 cm
15000 spinner bottles	16个	40个	80个	160个	400个	800个	1600个
T225 square bottles	160个	400个	800个	1600个	4000个	8000个	16000个
Single layer cell factory	56个	141个	283↑	566个	1415个	2830个	5660个

SPECIFICATION

Туре	Suspension cell culture bag	Perfusion membrane cell culture bag	Fabric Carrier Cell Culture Bag				
Total volume	2L/5L/10L/20L 50L/100L/200L	2L/5L/10L/20L 50L/100L/200L	2L/5L/10L/20L 50L/100L/200L				
Working volume	10%-50% of total volume	10%-50% of total volume	50% of total volume				
Carrier quantity/ surface area	No	No	30g/L, the specific surface of the carrier is 1200cm2/g				
Perfusion membrane	No	Yes	No				
Package	Double separate packing						
Sterilization	R-ray irradiation sterilization (25-40kGy)						
Temperature	-45 to +45 degrees						





SKC900 STIRRED BIOREACTOR



SKC900 single-use stirred bioreactor is including the excellent mixing and stirring effects of traditional stirred bioreactor, besides it's including single-use advantages of easy operation, no need for cleaning validation, and no need for sterilization, therefore It is an ideal equipment for production of vaccine, recombination protein and monoclonal antibody.

SKC900 single-use stirred bioreactor composed of a stainless steel container, a single-use cell culture bag and a controller equipped with excellent sensors.

SKC900 single-use stirred bioreactor can match the bags volume from 50L to 2000L, and is easy to achieve suspension cell culture and micro-carrier/fabric carrier adherent cell culture.

SKC900 single-use stirred bioreactor can simplify the scale-up process, and be easily converted between stainless steel and single use bioreactors, reduce the risk of bioprocess transfer.

FEATURE



Sterility/enclosed/low cost: Enclosed environment with low contamination risks and less maintenance



Convenient validation and operation: No need for cleaning and sterility validation with low contamination risks



Own culture bags: Cell culture bags are produced by ourselves with validation and stable long-term supply



Scale-up: 50L to 2000L available



Validation support: providing validation services of single-use cell culture bags and bioreactors



Monitoring: Monitoring cell culture process



Customization: Replace abroad bags or bioreactor bags



GMP standards: 3 level security/electronic signature/audit trail/GMP standards

APPLICATION

Suspension cell culture

Adherent cell culture on micro-carriers or fabric carriers

Culture of mammalian cells/insect/stem cells

Medium and low density microbial cell culture

Pilot scale production

R&D scale-up process

Large-scale seed culture

Production of vaccines, recombination proteins and antibodies

SPECIFICATION

Container size	50L	200L	500L	1000L	2000L			
Bag volume	50L	200L	500L	1000L	2000L			
Controller system	Siemens PLC, 15-inch Siemens color touch screen, Ethernet interface, USB interface.							
Temperature control	Water jacket, PID automatic control, Pt100 sensor, accuracy \pm 0.5°C, PID control, temperature range 8-45°C, Max. 60°C							
Max. speed	360RPM	350RPM	250RPM	140RPM	115RPM			
Gas*	4 thermal mass flow controller(Air/O2/N2/CO2), 1 thermal mass flow controller with 3 gas mixer and overlay optional							
Exhaust **	filter heater/0.2μm air filter							
Exhaust analysis	O2 (0-50%), CO2 (0-15%)							
Gas pressure control	Measured and controlled							
рН※	1.Single-use fiber pH probe, control range: 6-8, accuracy: ±0.05pH; PID control; cascaded by peristaltic pump/carbon dioxide concentration etc.							
	2.Repeatable use of traditional pH probe with sterilizable electrode adapter, control range: 4-10, PID control, cascaded by peristaltic pump, carbon dioxide concentration, etc.							
DO*	1. Single-use fiber DO electrode, range: 0-100%, PID control; cascaded by rocker speed/MFC etc.							
	2.Repeatable use of traditional DO probe with sterilizable electrode adapter, range: 0-100%,PID control, cascaded by rocker speed/MFC etc.							
Pumps*	4 fixed-speed pumps for addition/acid/base/anti-foam etc., 25PRM/100PRM variable-speed pumps optional							
Load cell ※	Optional							
Signal interface*	Optional							
Utilities	Air/O2/N2/CO2; Requirement:30-35 psig							
Power supply	AC220V,50/60Hz,20A							

Note: *Optional





PROCESS APPLICATION

PROCESS APPLICATION

VECTOR VIRUS PROCESS FLOW








Medium preparation


Medium balance

Inoculation seed cells









Circulating or perfusion cell culture

Cell growth monitoring

Inoculation virus or transfection

Batch/continuous harvesting virus

PERFUSION MEMBRANE N-1 SCALE-UP PROCESS FLOW









2x5mL centrifuge tubing (5*10'cells/ml)

Perfusion membrane cell culture (1L, 4*10°cells/ml)

Perfusion membrane perfusion cell culture (10L, 8*10' cells/ml)

2000L single-use bioreactor.

STEM CELL CULTURE PROCESS FLOW







Medium preparation

Medium balance

Inoculation cells









The supernatant

Cells

Circulating perfusion cell culture

Growth monitoring

Harvesting the supernatant/cells

IMMUNE CELL CULTURE PROCESS FLOW







Medium preparation

Medium balance

Inoculation cells







Perfusion membrane perfusion cell culture

Growth monitoring

Harvesting immune cells



VACCINE/VIRUS/CELL PRODUCTION SCALE-UP PROCESS FLOW

Scale-up process from fabric carrier cell culture to micro-carrier cell culture







Scale-up culture of 5L fabric carrier bag, trypsinization of adherent cells



14L/40L micro-carrier culture bioreactor

Scale-up process of fabric carrier bag cell culture and trypsinization



The fabric carrier shake flask culture, harvest the total number of (5-15)*10°cells, and inoculate it on a 5L fabric carrier bag, the inoculation density is (2-6)*10°cells/ml



5L culture bag grows for 5-7 days and digest, then inoculate it into 1-3 per 50L fabric carrier bags



50L culture bag grows for 5-7 days and digest, then harvests (1-1.5)*10¹¹ cells, therefore can inoculate multiple 200L fabric carrier bags or bioreactors, inoculation density is 6*10⁵ cells/ml



Multiple 200L fabric carrier bag/ bioreactor culture for production

293T CULTURE PROCESS OF FABRICE CARRIER CULTURE BAG




Fabric carrier shake flask 293T culture, and trypsinizated and then inoculate 3*10'cells



2L Culture bag grows for 5-7 days and trypsinizated, then take out 2*10¹⁰ cells inoculate into 50L culture bag



Culture bag 50L grow for 72 hours, change the medium, transfect, then change the medium for 6 hours, and culture for 90 hours by circulate perfusion



Harvest, Trypsinizated from fabric carrier and then lysed, or lysed by lysate





SKC300 BIOREACTOR



INTRODUCTION

Adopts brand components from domestic and abroad. The equipment has a neat appearance, easy to clean and sterilize, concise human-machine interface, stable process control and convenient operation.

It is a highly versatile bioreactor, compatible with most mainstream cell culture stirring blades and bioreactors.

FEATURE



Multiple Impellers

- Pitched blade impeller or spin filter type Impeller for suspension cell culture
- Packed-bed basket Impeller for adherent cell perfusion culture
- · Cell lift Impeller for Micro-carrier culture



Multiple water-jacketed vessel: 2.5L, 5.0L, 7.5L and 14.0L optional



Preset control algorithm: Preset pH, DO, gassing (four gas mixture and gas flow), temperature and others, easier to culture



Multiple cascades: Customized parameter cascade/time cascade/alarm correlation algorithms



Siemens PLC control system: Controlling stirring speed/ temperature/PH/DO/gassing/anti-foam/liquid level etc.



Modular design: Flexible and easy to maintain



Power-off self-recovery: Automatically recover the state before power-off



Industrial Ethernet: Control multiple devices remotely and portable devices optional

HMI



SIEMENS 12 INCH FULL-COLOR TOUCH SCREEN

Siemens 12-inch full-color touch screen is adopted to record and save process data in real time so as to meet GMP standards on drug production.



PARAMETER SETTINGS

Clear parameter setting screen, trend in real time and history, alarm setting, simple and friendly human-computer interface to show control process.




3 level security, electronic signature, audit trail, batch report, log record and GMP standards









Packed-bed basket Impeller for fabric carrier cells

PROCESS CHARACTERISTICS

- Larger adhesion area and strong cell adhesion ability
- No contact aeration system, low (no) shear forces stirring system
- No need to pre-treat the fabric carrier, and fast adherent
- Unrestricted perfusion culture without interception device
- Low cell seeding density, 3D fabric carrier supports ultra-high cell density culture
- · High viral titer, better for purification

PROCESS APPLICATION

- R&D and production for adherent CHO cells of EPO/TPO
- R&D production for human virus vaccines, e.g. COVID-19 vaccine/rabies vaccine/EV71 vaccine/IPV vaccine /influenza vaccine/vaccinia, etc.
- R&D and production for LV (lentivirus)/ ADV(adenovirus)/ AAV(adeno-associated virus)/ RV(retrovirus)/ HSV (oncolytic virus)/ poxvirus and other viral vectors of 293/ 293t culture
- R&D and production for influenza vaccine of MDCK
- R&D and production for veterinary virus vaccines, e.g. swine fever/PRRS/ swine diarrhea/ swine gastroenteritis/ avian influenza/ pseudorabies etc.
- R&D and production for cell exosomes
- · High-density culture of hybridoma cells/stem cells etc.
- R&D and production for primary cells and virus vaccines, e.g. human rabies vaccines from chicken embryos, etc.

Cell lift Impeller for Micro-carrier cells


PROCESS CHARACTERISTICS

- Larger adhesion area and higher density
- Scale up and cultivated
- Observe and count easier for sampling
- No contact aeration system, low (no) shear forces stirring system
- Better for perfusion cell culture
- Trypsinization harvest or scale-up culture

PROCESS APPLICATION

- R&D and production of EV71 (Enterovirus 71) vaccine for VERO
- R&D and production of human rabies vaccine for VERO
- Culture for stem cells, e.g. IPS and MSC etc.
- R&D and production for human rabies vaccines of diploid cells

PRODUCT MODEL



APPLICATION





03 Micro-carrier adherent cell culture for adherent cells





SPECIFICATION


Vessel volume	Total volume	2.5L	5.0L	7.5L		14.0L	
	Working volume	0.75L-1.75L	1.25-3.75L	2.0-5.5L		3.0L-10.5L	
Headplate connector	Qty.*Size	1*6.35mm 9*PG13.5	3*6.35mm、1*19mm 10*PG13.5	3*6.35mm、1* 12*PG13.5	19mm	3*6.35mm、1*19mm 12*PG13.5	
Control system	Hardware *	Siemens PLC, 12-inch Siemens touch screen, 4 pumps, Ethernet interface, US interface, load cell optional					
Temperature control	NUMBER OF THE PROPERTY OF THE		neater PID automatic control, Pt100 sensor, accuracy ±0.1°C, temperature B-45°C above coolant temperature, Max.60°C				
Agitating system	Agitating form	Upper agitating, magnetic-driven					
	Impeller*	Pitched-blade impeller/basket impeller/cell lift impeller/spin filter impeller optional					
	Agitating speed	Speed range:25-500RPM, accuracy: ±1RPM (excluding 14.0L) Speed range			ge: 25-150RPM (14.0L)		
Exhaust	Condensate/ filter	Stainless steel/0.2µm air filter					
Gas ※	Control system	1 or 4 thermal mass flow controllers optional					
Sparge ※	Circinate. Mic	cro and macro sparge optional					
рН፠	2-12, PID cont	D control cascaded with pumps/gas flow. ISM probe optional					
DO%	0-200%, PID c	-200%, PID control cascaded by agitation/pumps/gas flow. ISM probe/ optical probe optiona					
Pumps*	Parameters	4 fixed-speed pumps for addition/acid/base/anti-foam etc., speed 25PRM, 25PRM /100PRM fixed-speed or variable-speed pumps optional.					
Utilities	Water/gas	Air/O2/N2/CO2, 10PSIG					
Power supply	AC220V, 50/60Hz, 10A						

Note: *Optional





SKC800 BIOREACTOR





N/A connector



Sampling valve



Temperature validation

INTRODUCTION

A bioreactor that can automatically sterilize.

Adopts self-developed biological reaction process control algorithm. By stirring, temperature, pH, DO, sparger, overlay, peristaltic pump, anti-foam and liquid level, load cell, cold spot detection and additional analogy input and output to achieve batch/fed batch/continuous process.

Adopts brand components from domestic and abroad. The equipment has a neat appearance, easy to clean and sterilize, concise human-machine interface, stable process control and convenient operation.

A bioreactor that compatible with cell culture stirring Impeller and bioreactors, e.g. Suspension cell culture suits pitched blade impeller or spin filter type Impeller; Perfusion culture of fabric carrier adherent cells suits packed-bed basket Impeller; Micro-carrier cell perfusion culture suits cell lift Impeller.

HARDWARE



Well-known brand components, e.g. 4-valve/sampling valve/ bottom valve/mass flow controller, etc.



N/A connector and ferrule fitting




Cold point monitoring and F0 of sterilization to meet GMP standards




Multiple stirring blades available, e.g. batch cell culture/fed-batch cell culture/continuous perfusion culture/cell suspension culture/high-density culture of micro-carriers and fabric carriers

HMI

- · Siemens PLC: Running stability
- Siemens color touch screen: To record and save process data and human-machine interface simple





- Multiple cascades: Customized, e.g. parameter cascade/time cascade/alarm correlation algorithms
- Industrial Ethernet: Control multiple devices remotely and portable devices optional
- Modular design: Flexible and easy to maintain.
- 3 level security, electronic signature, audit trail, batch report, log record and GMP standards




4 valve group



Cold spot monitoring



Packed-bed basket Impeller with fabric carrier



Flow chart

Simple process interface, including system activation/four-gas control/temperature/stirring/peris -taltic pump/pH DO control, etc.



Automatic Sterilization Interface —

Directly set, e.g. drain time/sterilization time/sterilization temperature/culture temperature etc., and the on-off status of control valve can be monitored



Audit

Audit trail, batch report, log record and GMP standards





Packed-bed basket impeller for fabric carrier cells

PROCESS CHARACTERISTICS

- No contact aeration system, low (no) shear forces stirring system
- No need to pre-treat the fabric carrier, and fast adherent
- Unrestricted perfusion culture without interception device
- Low cell seeding density, 3D fabric carrier supports ultra-high cell density culture
- High viral titer, better for puritification

PROCESS APPLICATION

- Larger adhesion area and strong cell adhesion
 R&D and production for adherent CHO cells of EPO/TPO
 - R&D and production for human virus vaccines, e.g. COVID-19 vaccine/rabies vaccine/EV71 vaccine /IPV vaccine /influenza vaccine/vaccinia, etc.
 - R&D and production for LV (lentivirus)/ADV (adenovirus)/AAV(adeno-associated virus)/ RV (retrovirus) /HSV (oncoly -tic virus) /poxvirus and other viral vectors of 293/293t cell
 - R&D and production for influenza vaccine of MDCK
 - R&D and production for veterinary virus vaccines, e.g. swine fever/PRRS/swine diarrhea /swine gastroen -teritis/avian influenza/pseudorabies etc.
 - R&D and production for cell exosomes
 - High-density culture of hybridoma cells/ stem cells etc.
 - R&D and production for primary cells and virus vaccines, e.g. human rabies vaccines from chicken embryos, etc.



Cell lift impeller for micro-carrier cells

PROCESS CHARACTERISTICS

- Larger adhesion area and higher density
- Scale up and cultivated
- Observe and count easier for sampling
- No contact aeration system, low (no) shear forces stirring system
- Better for perfusion culture
- Trypsinization harvest or scale-up culture

PROCESS APPLICATION

- R&D and production of EV71 (Enterovirus) 71) vaccine for VERO
- R&D and production of human rabies vaccine for VERO
- Culture for stem cells, e.g. IPS and MSC etc.
- R&D and production for human rabies vaccines of diploid cells

Pitched blade impeller for suspension and low-density micro-carrier cells

PROCESS CHARACTERISTICS

- Easy to scale up for suspension and low -density micro-carrier cell culture
- Observe and count easier for sampling
- Flexible process for batch/fed-batch/continuous perfusion cell culture
- Tolerate shear force of stirring Impeller and I nfluence of air bubbles
- Together with ATF/spin filter, etc. to achieve continuous perfusion cell culture

PROCESS APPLICATION

- Protein of CHO cell antibodies
- Suspended animal cells culture
- Human rabies vaccine of VERO cells



APPLICATION



For fabric carrier adherent culture, especially for secreted products, e.g. CHO for EPO/TPO, VERO cells for human rabies vaccine/Japanese encephalitis vaccine/COVID-19 vaccines/MDCK cells for flu vaccines, etc.



02 For micro-carrier adherent culture, e.g. VERO cells/diploid cells for human rabies vaccine, etc.



03 For suspension cell culture, e.g. CHO cell antibody/SF9 cell cervical cancer vaccine, etc.

SPECIFICATION

Vessel volume	Vessel volume	40L	75L	150L	300L			
	Working volume	10.75L-32.0L	18.75-60.0L	37.5-120.0L	75.5L-240.0L			
Control system	Siemens PLC, Siemens touch screen, 3 pumps, Ethernet interface, USB interface							
Temperature control	Water jacket, electrical-heater PID automatic control, Pt100 sensor, control accuracy ±0.1℃, PID control, temperature control range 8-45℃ above, Max. 60℃							
Agitating system *	Agitating form	n Upper mixing, double seal						
	Impeller*	Pitched-blade impeller/basket impeller/cell lift impeller/spin filter impeller optional						
	Agitating speed	25-200rpm	25-150rpm	25-120rpm	25-100rpm			
Load cell	Load cell mould	0-40KG, cascaded by peristaltic pump optional	0-100KG, cascaded by peristaltic pump optional	0-200KG, cascaded by peristaltic pump optional	0-300KG, cascaded by peristaltic pump optional			
Sparger ※	4 gas mixing, 1 or 4 thermal mass flow controllers optional							
Overlay ※	4 gas mixing, 1 or 4 thermal mass flow controllers optional							
Exhaust	Exhaust condenser, exhaust filter heater							
Gas distributor ※	Circinate. Micro and macro sparge optional							
рНЖ	2-12, PID control cascaded by pumps/gas flow. ISM probe optional							
DO*	0-200%, PID control, cascaded by agitation/pumps/gas flow. ISM probe/optical probe options							
Sterilization	Automatic SIP							
Cold spot detection ※	Optional							
Pumps ※	Fixed-speed pumps for addition/acid/base/anti-foam etc., speed 100PRM,							
Level and anti- foam ※	1 level electrode. High level electrode optional							
SCADA software **	Optional. GMP/GAMP standards.							
Utilities	Water: 50um filtration; industrial steam 30-35PSIG, pure steam 30-35PSIG; 4 gases 30-35PSIG							
Power supply	AC220V,50/60HZ,20A							

Note: * Optional



PRODUCT INTRODUCTION

SKC700 BIOREACTOR

INTRODUCTION

A cost-effective culture system. It consists of a master controller and quality sterilization glass vessel.

By controlling the agitation/temperature/pH/ DO/sparger/pumps/ anti-foam and level probe to achieve plant and animal cell culture, microbial fermentation.

Adopts brand components from domestic and abroad. The equipment has a neat appearance, easy to clean and disinfect, concise human-machine interface, stable process control and convenient operation.



Parallel Control Software

Parallel control and record 2 to 8 bioreactors parameter, 3 level security, audit trail, batch report, log record and GMP standards.

FEATURE

- Multiple water-jacketed vessel: 2.5L, 5.0L, 7.5L and 14.0L optional
- Multiple cascades: Customized parameter cascade/time cascade/alarm correlation algorithms
- Siemens PLC control system: For stirring speed/temperature/pH/DO gassing/anti-foam/liquid level etc.
- A 7-inch color industrial-grade touch screen: To record and save process data
- Industrial Ethernet: Control and record 2 to 8 bioreactors remotely
- 3 level security, electronic signature, audit trail, batch report, log record and GMP standards
- Modular design: Flexible and easy to maintain
- Power-off self-recovery: Automatically recover the state before power -off



Note: *Optional

APPLICATION

- R&D and cultivation for antibody protein drugs
- R&D and cultivation for stem cells
- R&D and cultivation for human or veterinary vaccines
- R&D and cultivation for immune cells

CELL CULTURE BIOREACTOR CONFIGURATION



PRODUCT MODEL

Item	Model	Heating Type	Applicable Scope
SKC700 Bioreactor Fermenter	SKC0701-1401	Blanket	Applicable to SKC0700 series, 2.5L vessel kit. Magnetic drive
	SKC0701-1402	Blanket	Applicable to SKC0700 series, 5L vessel kit. Magnetic drive
	SKC0701-1405	Blanket	Applicable to SKC0700 series, 7.5L vessel kit. Magnetic drive
	SKC0701-1410	Blanket	Applicable to SKC0700 series, 14.0L vessel kit. Magnetic drive

SPECIFICATION

	Vessel volume	2.5L	5.0L	7.5L	14.0L				
Vessel volume	Working volume	0.75-1.75L	1.25-3.75L	2.0-5.5L	3.0-10.5L				
	Fermenter type	Heating blanket	Heating blanket single-wall, borosilicate glass						
Headplate connector	Qty.*Size	1*6.35mm 9*PG13.5mm	3*6.35mm 10*PG13.5mm 1*19mm	3*6.35mm 12*PG13.5mm 1*19mm	3*6.35mm 12*PG13.5mm 1*19mm				
Control system	Hardware **		Siemens PLC, 7-inch color touch screen, 3 pre-pumps, Ethernet interface, USB interface, load cell optional						
Temperature control *		tet temperature control,PID automatic control, Pt100 sensor, accuracy ±0.1°C, control range 8-45°C, Max.60°C							
2 2	Agitating form	Upper agitating, magnetic-driven							
Agitating system	Impeller*	Pitched-blade impeller/basket impeller/cell lift impeller/spin filter impeller optional							
2,000	Agitating speed	Speed range: 25-500RPM,accuracy: ±1RPM (excluding 14.0L) Speed range:25-150RPM (14.0L)							
Exhaust	Condensate/filter	Stainless steel/0.2µm air filter							
Canif	Control system *	1 or 4 thermal mass flow controllers optional							
Gas*	Sparge **	Circinate. Micro and macro sparge optional							
рНЖ	2-12, PID control cascaded by pumps/gas flow. ISM probe optional								
DO%	0-200%, PID control, cascaded by agitation/pumps/gas flow. ISM probe/optical probe optional								
Pumps*	3 fixed-speed pu pumps optional	B fixed-speed pumps for addition/acid/base/anti-foam,etc.,speed 25PRM, variable speed pumps optional							
Utilities	Water/gas	Air/O2/N2/CO2, 10PSIG							
Power supply	AC220V, 50/60Hz, 10A								



APPLICATION

- Plasmids, human proteins, etc. in Escherichia coli
- Pichia pastoris to induce the production of microbial fermentation proteins, etc.
- · Fermentation culture of streptomycin yeast/baculovirus
- Batch/fed perfusion culture of Escherichia coli, yeast or plant cells



MICROBIAL FERMENTOR CONFIGURATION



Note: *Optional

MICROBIAL FERMENTER MODEL

Item	Model	Heating Type	Applicable Scope
SKC700 Microbial Fermente	SKC0701-1151	Blanket	Applicable to SKC0700 series, 2.5L vessel kit. Direct drive
	SKC0701-1152	Blanket	Applicable to SKC0700 series, 5L vessel kit. Direct drive
	SKC0701-1155	Blanket	Applicable to SKC0700 series, 7.5L vessel kit. Direct drive
	SKC0701-1150	Blanket	Applicable to SKC0700 series, 14.0L vessel kit. Direct drive

SPECIFICATIONS

	Vessel volume	2.5L	5.0L	7.5L	14.0L			
Vessel volume	Working volume	0.75-1.75L	1.25-3.75L	2.0-5.5L	3.0-10.5L			
	Fermenter type	Heating blanke	Heating blanket single-wall, borosilicate glass					
Headplate connector	Qty.*Size	1*6.35mm 9*PG13.5mm	3*6.35mm 10*PG13.5mm 1*19mm	3*6.35mm 12*PG13.5mm 1*19mm	3*6.35mm 12*PG13.5mm 1*19mm			
Control system	Hardware*		Siemens PLC, 7-inch color touch screen, 3 pre-pumps, Ethernet interface, USB interface, load cell optional					
Temperature control *		et temperature control,PID automatic control, Pt100 sensor, accuracy ±0.1°C, ontrol range 8-45°C, Max.60°C						
	Agitating form	Upper agitating, magnetic- driven						
Agitating system	Impeller*	Rushton blade						
System	Agitating speed	Speed range: 25-1200RPM, accuracy: +1R PM						
Exhaust	Condensate/filter	Stainless steel/0.2µm air filter						
C\(\forall \)	Control system *	1 or 4 thermal mass flow controllers optional						
Gas**	Sparge *	Circinate. Micro and macro sparge optional						
рН፠	2-12, PID contro	2-12, PID control cascaded by pumps/gas flow. ISM probe optional						
DO%	0-200%, PID cont	-200%, PID control, cascaded by agitation/pumps/gas flow. ISM probe/optical probe optional						
Pumps*		d-speed pumps for addition/acid/base/anti-foam,etc.,speed 25PRM, 25PRM variable- d pumps optional						
Utilities	Water/gas	Air/O2/N2/CO2, 10PSIG						
Power supply	AC220V, 50/60Hz, 10A							

SERVICES

Wuhan CEKG Technology CO., LTD. provides customers with best on-site services based on its service mission.

SERVICES

Wuhan CEKG Technology Co., Ltd. has own cell culture laboratory and technical support team, which can provide cell culture technical services, e.g. process scale-up and validation, process scale-up, process optimization, etc., and can provide a whole cell culture solutions.





SERVICES

Wuhan CEKG Technology Co., Ltd. provides installation and test, technical consultation, maintenance, and IQ/OQ/PQ validation services for different bioreactors and liquid distribution and storage systems.

Single-use product validation services

- · Chemical compatibility · Extractables
- Leachable
- Adsorption



SERVICES

Wuhan CEKG Technology Co., Ltd. has various spare parts with bioreactors and consumables with liquid mixing/storage system.

Spare parts: Peristaltic pumps/solenoid valves/diaphragm valves/water pumps/mass flow controllers/Siemens PLC module/pH or DO probes/DO probe membranes or electrolytes/glass joints/single-wall or double-wall glass vessels/sampling kits/O-rings/gaskets/diaphragms/various Impeller kits/316L meshes, etc.

Consumables: CPC quick connectors/CPC aseptic connectors/T-shape or Y-shape or 4 ways connectors/reducing T connectors/straight connectors/silicone tubing/C-flex tubing/single-use addition bottle/silicone tubing clamp, etc.





△ pH and DO probes



△ Peristaltic pump



△ Connector accessories



△ O-rings, gasket, diaphragm



△ CPC connectors




△ Circulating water pump



△ Silicone tubing and C-flex tubing





Wuhan CEKG Technology Co., Ltd.

- ♦ ADD: No. 219, Guanggu new power, Guanggu Second Road, Donghu New Technology Development Zone, Wuhan
- **28** 86 (027) 87688776 86 (027) 65521996
- www.whcekg.com
- whcekgsales6@163.com whcekg@163.com

